



Made in Germany

# Suflex™

Hydrophilic VPS Impression Material

**HIOSSSEN**

- ✓ Excellent thixotropic properties
- ✓ Clear Visiability
- ✓ Hydrophilic properties
- ✓ Cost Effective



# Suflex™ Hydrophilic VPS Impression Material



Suflex™ is an A-Silicone based impression material that has remarkable elastic properties and high tensile strength. The combination of low viscosity allows for optimal flow and high thixotropy (non-drip, no sagging) make manipulating Suflex easy and convenient. The Suflex™ hydrophilic and wettable properties assist with forming void-free impressions with superior margin reproduction and capture accurate details even in the presence of saliva.



**VPS (Vinyl Polysiloxane)** : Modern Impression material based on A-Silicone

## Suflex™ Features & Benefits

Features	Benefits
<ul style="list-style-type: none"> <li>• Thixotropic</li> <li>• Clear Visiability</li> <li>• Hydrophilic</li> </ul>	<ul style="list-style-type: none"> <li>- Easy and convenient manipulation (no dripping or sagging, stays in the tray)</li> <li>- Suflex offers the clear visiability by presenting the blue-colored margin on the yellow field after impression.</li> <li>- It's affinity to water, make it wettable helps create void-free and accurate impressions.</li> </ul>

## Suflex™ Mixing Time

	Light 50ml	Mono 50ml	Heavy 50ml	Putty 900ml	Auto Mono 362ml	Auto Heavy 362ml
<b>Working Time (min.)</b>	≤ 2:00	≤ 1:30	≤ 2:00	≤ 2:00	≤ 2:00	≤ 2:00
<b>Setting Time in the mouth (min.)</b>	≥ 2:00	≥ 2:30	≥ 2:00	≥ 2:00	≥ 2:00	≥ 2:00
<b>Color</b>						

## Suflex™ Order Information

PRODUCT	CODE	Quantity	COLOR
Suflex Putty	SUFPA	450ml * 2	
Suflex Light	SUF50LA	50ml * 2	
Suflex Light Econo refill	SUF50LB	50ml * 10	
Suflex Mono	SUF50MA	50ml * 2	
Suflex Mono Econo refill	SUF50MB	50ml * 10	
Suflex Heavy	SUF50HA	50ml * 2	
Suflex Heavy Econo refill	SUF50HB	50ml * 10	
Suflex Auto Heavy Intro Pack	SUFAHS	362ml * 1	
Suflex Auto Heavy Econo refill	SUFAHA	362ml * 2	
Suflex Auto Mono Intro pack	SUFAMS	362ml * 1	
Suflex Auto Mono Econo refill	SUFAMA	362ml * 2	

